

1

Assemble Sizing Chart First


2011

Important

Before printing, turn off "fit to page" in your print options dialogue box.*

A

← 1½" →



B

Fold under on *Solid Fold Line* on page 2, align *Solid Lines* and arrows, and tape pages together.

2


Length Measurement

A

- Fold sizing chart on dotted line
- Make a 90-degree angle
- Place on hard surface
- Put fold firmly against wall
- Tape to prevent mis-sizing

use arrows to align pages

Solid Line



children to adult sizes

use arrows to align pages

use
arrows
to align
pages

B

- Place heel of AFO on sizing chart
- Be certain to include posterior stop
- If foot is longer than brace, include foot

C

- Trace AFOs for accuracy
- Match longest point to number on right
- If you are between or slightly above line, round up


D IMPORTANT

Necessary Adjustments for Circumference

For all AFO's, TCFlex, DAFO's, KAFO's, SMO's with:

- Solid ankle, instep strap, articulating hinges, posterior stop
Increase 1 full size
- Turbo boot, posterior stop, hinges, footplates, instep, forefoot and toe strapping:
Increase 2 sizes

Fold under on *Solid Fold Line*,
align *Solid Lines* and arrows,
and tape pages together.



Please be sure to accurately measure AFOs using this sizing chart and follow steps A–D listed on left. For best results, please verify square in left hand corner measures 1½" wide x 1½" high. KeepingPace cannot guarantee 100% accuracy due to discrepancies in various printing capabilities.

Fold up on *Dotted Fold Line*.



KEEPINGPACE
Children's Orthopedic Footwear

Keeping Pace, Inc.

contact us:
footwear@keepingpace.com
P (888) 526-0020
F (978) 526 0025
USA

visit us on the web:
www.keepingpace.com

use
arrows
to align
pages