


2011

**Important**

Before printing, turn off "fit to page" in your print options dialogue box.\*



**Length Measurement**

- 1**
  - Place bottom edge of this sheet against wall on a hard surface.
  - Tape down to prevent any mis-sizing
- 2**
  - Place heel of AFO on sizing chart **against wall**
  - Be certain to include posterior stop
  - If foot is longer than brace, include foot
- 3**
  - Trace AFOs for accuracy
  - Match longest point to number on right
  - If you are between or slightly above line, round up

**4 IMPORTANT**

**Necessary Adjustments for Circumference**

For all AFO's, TCFlex, DAFO's, KAFO's, SMO's with:

- Solid ankle, instep strap, articulating hinges, posterior stop  
*Increase 1 full size*
- Turbo boot, posterior stop, hinges, footplates, instep, forefoot and toe strapping:  
*Increase 2 sizes*



**children and adolescent sizes**



Keeping Pace, Inc.  
*contact us:*  
 footwear@keepingpace.com  
 P (888) 526-0020  
 F (978) 526 0025  
 USA  
*visit us on the web:*  
 www.keepingpace.com